

**METHODS OF DETERMINING DEFICIENCIES IN INTRACELLULAR
LEVELS OF CYSTEINE AND GLUTATHIONE**

BACKGROUND OF THE INVENTION

Cross-Reference to Related Applications

This is a divisional of U.S. Serial No. 10/017,625, filed
December 13, 2001, which is a divisional of U.S. Serial No. 08/922,279,
filed September 3, 1997, now abandoned.

Field of the Invention

The present invention relates generally to the fields of
nutrition and biochemistry and cellular glutathione metabolism. More
specifically, the present invention relates to measurement of levels of

intracellular function of cysteine and glutathione so as to provide one measurement an individual's capability of preventing degenerative disease and dealing with oxidative stress.

5 Description of the Related Art

It is now accepted widely that a number of human health conditions, including aging, arthritis, cancer, atherosclerosis, myocardial infarction, stroke, viral infection, pulmonary conditions, bowel disease and neurodegenerative disease, can develop or be worsened by the
10 presence of reactive oxygen molecules, commonly referred to as free radicals. These hostile molecules are normal by-products of physiological processes and are produced by metabolism of oxygen; e.g., via cellular respiration or immune system function (killing of foreign materials), and by numerous enzymatic reactions essential for
15 metabolism. In addition, free radicals are found commonly in the environment. Environmental sources of free radicals include smoke, ionizing radiation, air pollution, chemicals (carcinogens, many petrochemicals, biocides, dyes, solvents, cytostatic drugs, etc.), toxic heavy metals and oxidized or rancid fats. Some of the most common
20 free radicals are super oxide, hydroxyl radical, singlet oxygen, and peroxides, including hydrogen peroxides. Certain valences of iron and

copper can catalyze formation of free radicals, which, although short-lived, promote a chain reaction of radical formation, followed by a wake of altered, damaged biological molecules.

Free radicals are toxic to living organisms, causing structural damage to biological molecules. Molecular damage may result in alteration of genetic codes, disruption of cell membrane integrity, neurological disorders, endocrine imbalances, increased allergies, vascular endothelial destruction and joint degradation and inflammation.

Protection from the deleterious effects of free radicals is found in a diverse range of molecules termed antioxidants. Free radicals and their chain by-products can be neutralized and converted to less harmful products by antioxidants. Antioxidants may be enzymes (such as superoxide dismutase, catalase, glutathione peroxidase), essential nutrients (such as beta carotene, vitamins C and E, selenium) or a wide variety of endogenous compounds (such as glutathione) or dietary compounds (such as the bioflavanoids). Thus, the human body has several, natural quenchers of free radicals.

Research in humans has indicated that deficient intakes of nutrient antioxidants are associated with higher risks of cancer, cardiovascular disease, arthritis, cataracts, etc. Also, a higher intake of

nutrient antioxidants is associated with lower incidence of chronic degenerative diseases. Encouraging studies indicate that intervention with antioxidant nutrient supplements may have therapeutic benefit in humans.

5 Laboratory analysis of antioxidant status has not become routine for a variety of reasons. Free radicals are extremely fleeting and generally are not amenable to direct measurement. By-products of free radical damage can be measured, such as malondialdehyde (MDA), thiobarbituric acid-reactive substances (TBARS) or lipid peroxides in
10 serum or urine; however, these tests may be indicators of oxidative stress but reflect only damage to certain types of biomolecules (mostly polyunsaturated lipids and nucleic acids). Yet methods of measurement of antioxidant nutrient levels in serum or cells and activities of antioxidant enzymes in cells could identify deficient levels of specific
15 components.

 Glutathione is a preeminent cellular antioxidant substance which is abundant in cytoplasm, nuclei and in mitochondria. In addition to glutathione's powerful antioxidant function, it serves also as a powerful anti-toxin, allowing the body to eliminate numerous xenobiotic
20 and carcinogenic compounds. Further, glutathione is essential for cell-mediated immune functions and it is critical to the maintenance of the

integrity of red blood cells. Moreover, it is recognized generally that deficiencies in the glutathione system lead to significant cellular aging, and, ultimately, cellular morbidity.

The concentration of cellular glutathione has a significant
5 effect on antioxidant function; and nutrient limitation, exercise and oxidative stress have significant effects on cellular glutathione concentrations. Under oxidative conditions, glutathione function can be depleted considerably through conjugation to xenobiotics, and by secretion of glutathione conjugates and glutathione disulfide from the
10 affected cells. A considerable amount of glutathione may become protein bound during severe oxidative stress. Fortunately, however, compounds such as N-acetyl-cysteine are available to increase intracellular glutathione function.

Glutathione is synthesized in a series of biochemical
15 reactions utilizing ATP, magnesium and the three amino acids glycine, glutamate and cysteine. In general, the rate of synthesis of gamma-glutamylcysteine determines the rate of synthesis of glutathione, and the sulfhydryl group of cysteine provides glutathione with its biological potency. Thus, measurement of cysteine availability is essential in
20 determining the functional availability of glutathione and in the assessment of antioxidant function.

Assessment of cysteine and glutathione also is helpful in assessing selenium deficiency. When glutathione functions as an antioxidant, it reacts with hydrogen peroxide to form glutathione disulfide in a reaction catalyzed by glutathione peroxidase. Glutathione peroxidase requires selenium as a functional cofactor. Thus, adequate cysteine function, combined with deficient or average SPECTROX™ results, is indicative of an intracellular selenium deficiency and calls for further testing.

The prior art is deficient in the lack of simple cost-effective means of assessing biochemically the levels of intracellular function of cysteine and glutathione in a human, and thus an individual's capability of dealing with oxidative stress. The present invention fulfills this longstanding need and desire in the art.

SUMMARY OF THE INVENTION

One object of the present invention is to provide a cell culture medium useful for determining a levels of intracellular function of glutathione in lymphocytes and for performing biochemical analysis of antioxidant function, said medium comprising a buffered, serum-free

solution containing the following ingredients: a carbohydrate selected from the group consisting of glucose or a compound biologically capable of producing glucose in the cells; a biologically usable form of pantothenic acid, choline or a biologically utilizable form of a substance
5 capable of producing choline in the cells; inorganic ions comprising chloride, phosphate, calcium, magnesium, potassium, sodium, and iron in a biologically utilizable form; L- Buthionine-[S.R.]-Sulfoximine; deionized water; and a mitogen in an amount effective to stimulate the lymphocytes being assayed; said buffered, serum-free solution having a
10 pH from about 6.8 to 7.6; and wherein said cell culture medium is characterized by being effective to determine an levels of intracellular function of glutathione in lymphocytes, nutritional deficiencies, inadequacies, and imbalances and to analyze biochemically antioxidant function of the lymphocytes.

15 One embodiment of the present invention includes supplementing the cell medium of the present invention with a nutrient supplement selected from the group consisting of biological utilizable forms of amino acids and vitamins, the nutrient being tested for being omitted from or being present in limiting or inhibitory amounts in the
20 nutrient supplement.

Another object of the present invention is to provide a method of determining a levels of intracellular function of glutathione and analyzing biochemically cellular antioxidant function in an individual, comprising the steps of: inoculating the cell culture medium
5 of the present invention with lymphocytes from an individual to be tested; incubating the inoculated cell culture medium; and comparing the response of the lymphocytes with an average response of lymphocytes from a control group of individuals.

A further object of the present invention is to provide a cell
10 culture medium useful for determining the level of intracellular function of cysteine and performing biochemical analysis of antioxidant function in human lymphocytes, said medium comprising a buffered, serum-free solution containing the following ingredients: a carbohydrate selected from the group consisting of glucose or a compound biologically capable
15 of producing glucose in the cells; a biologically utilizable form of pantothenic acid, choline or a biological usable form of a substance capable of producing choline in the cells; inorganic ions comprising chloride, phosphate, calcium, magnesium, potassium, sodium, and iron in a biologically utilizable form; cumene hydroperoxide; deionized
20 water; N-Acetyl-L-Cysteine; and a mitogen in an amount effective to stimulate the lymphocytes being assayed; said buffered, serum-free

solution having a pH from about 6.8 to 7.6, said cell culture medium characterized by being effective to determine nutritional deficiencies, inadequacies, and imbalances and to analyze biochemically antioxidant function of the lymphocytes.

5 One embodiment of the present invention includes supplementing the cell medium of the present invention with a nutrient supplement selected from the group consisting of biologically utilizable forms of amino acids and vitamins, the nutrient being tested for being omitted from or being present in limiting or inhibitory amounts in the
10 nutrient supplement.

 An additional object of the present invention is to provide a method of determining a level of intracellular function of cysteine and analyzing biochemically cellular antioxidant function in an individual comprising the steps of: inoculating the cell culture medium of the
15 present invention with lymphocytes from an individual; incubating the inoculated cell culture medium; and comparing the response of the lymphocytes with an average response of lymphocytes from a control group of individuals.

 Other and further aspects, features, and advantages of the
20 present invention will be apparent from the following description of the

presently preferred embodiments of the invention. These embodiments are given for the purpose of disclosure.

5

DETAILED DESCRIPTION OF THE INVENTION

The present invention is directed to a method useful for the biochemical analysis of levels of intracellular function of cysteine and glutathione in human lymphocytes. Such an analysis reflects the physiological health of an individual's antioxidant systems within the peripheral lymphocytes. The methodology of the present invention allows precise assessment of the intracellular function of cysteine and glutathione in human lymphocytes, so that therapeutic measures can be taken to improve an individual's antioxidant profile.

15

The present method, utilizing lymphocytes, offers distinct advantages, as lymphocytes: (1) are host to the cell-mediated immune system and are easily stimulated to grow (mitogenesis); (2) reflect time-averaged, long-term nutrient status (the life of a lymphocyte is about six months); (3) possess metabolic pathways common to other cells, and (4) are easily collected by standard venipuncture.

20

By measuring lymphocyte growth to assess intracellular function of cysteine and glutathione, the method of the present invention reflects the unique profile of each patient, which varies widely. Therapeutic treatment, therefore, can be tailored to the specific
5 biochemical requirements of the individual rather than the "average" patient as determined by so-called norms.

It will be apparent to one skilled in the art that various substitutions and modifications may be made to the invention disclosed herein without departing from the scope and spirit of the invention. The
10 following examples are given for the purpose of illustrating various embodiments of the invention and are not meant to limit the present invention in any fashion.

15 **EXAMPLE 1**

Drawing Patient's Blood

Two 10 ml whole blood specimens preserved with acid-citrate-dextrose are required for the method of the present invention.
20 No fasting is required. All that is required are blood drawing supplies. The assay of the present invention may be performed on the patient's

blood should be shipped at room temperature to a suitable laboratory. No centrifuging of the blood is required. Comprehensive test results provided a sound, scientifically-based analysis of the patient's glutathione or cysteine profile.

5

EXAMPLE 2

Sample Processing: Cell Isolation

10 All procedures are conducted using sterile techniques under a laminar flow hood to insure sample sterility. Each patient's blood samples are assigned an Accession number upon receipt at the laboratory. This accession number is used as the sample number to enable tracking throughout processing, data collection and data analysis
15 steps. Every test tube, centrifuge tube, microtiter plate and data printout involved with processing the patients sample is labeled with this sample (accession) number.

Each patient sample consists of (2) two Acid-Citrate-Dextrose (yellow top) vacutainer type tubes, each containing 8 ml of whole
20 blood. After being assigned an accession (sample) number, the whole

blood was mixed by inverting 6 times. The two tubes of whole blood were combined into a 50 ml disposable centrifuge tube.

A 500 μ l aliquot was removed aseptically from each sample and placed in a 12 x 17 mm tube. This aliquot was used to perform a whole blood cell count on the Coulter Cell Counter, Model T540. The whole blood cell count printout from the Coulter was labeled with the accession number and attached to the Worksheet for that patient.

Two (2) Ficoll gradient tubes were prepared for each sample by the addition of 5.0 ml of Histopaque 1077 (Ficoll/ Sodium Diatrizoate, Sigma Chemicals, St. Louis, Missouri) to each 15 ml conical centrifuge tube. Using a 10 ml pipette and an electric pipette aid, 8 ml of whole blood is slowly layered onto each of the Ficoll gradient tubes. The Ficoll gradient tubes were capped and centrifuged at 2160 RPM for 20 minutes.

After centrifugation was complete, the gradient tubes were carefully removed from the centrifuge to avoid disrupting the gradient. The buffy coat (containing the lymphocytes) found at the interface of the middle Ficoll layer was transferred using a 5 ml pipette into a 15 ml disposable conical centrifuge tube. The buffy coat was combined with phosphate buffered saline-0.72% glucose solution (PBS-G) to a final

volume of 12 ml. The tube was capped and inverted 6 times to mix the buffy coat and the PBS-G.

The tubes containing buffy coat and PBS-G were centrifuged at 2160 RPM for 5 minutes. After centrifugation, the supernatant was aspirated from the cell pellet and discarded. The cell pellet was resuspended into 12 ml of PBS-G, then inverted 6 times to insure adequate dispersal of the cell pellet. The sample was then centrifuged again as described above.

After the second centrifugation, the supernatant was aspirated and discarded. The cell pellet was resuspended in 6.0 ml of PBS-G. The cell pellet was disrupted and mixed with the PBS-G using a 5 ml pipette attached to an electric pipet aide. After a homogeneous cell suspension has been attained, a 200 μ l aliquot of the suspension was transferred into a 12 x 75 mm tube. This aliquot was used to perform the initial cell suspension (ICS) count with the Coulter Cell Counter Model T540.

The printout from this aliquot was labeled with the sample number and attached to the Worksheet. If the lymphocyte number was between 3.9 and 1.2 thousand cells per cubic millimeter (THSD/ mm^3) the sample is ready for plate inoculation. The volume of cell suspension to be added is found in TABLES I and IV. If the lymphocyte number was

greater than 3.9 THSD/ mm³, the sample must be rediluted. If, however, the lymphocyte number was less than 1.2 THSD/ mm³, the sample was rejected.

The amount of additional PBS-G that is required for proper
5 redilution was determined using the following calculations:

$$C_1 V_1 = C_2 V_2 \quad C = \text{Lymphocyte Concentration (THSD/mm}^3\text{)}$$

$$V = \text{Volume}$$

Where $C_2 = 3.0 \text{ THSD/mm}^3$, this was the desired lymphocyte
10 concentration of the final cell suspension. For example, when the initial
Cell Suspension Count from 6.0 ml = (LY# of 5.6 THSD/ mm³)

$$C_1, V_1 = C_2 V_2 \quad (5.6)(6.0 \text{ ml}) = (3.0)(X)$$

$$X = 11.2 \text{ ml final volume.}$$

15 An appropriate volume of PBS-G was added to resuspend cells
for a final volume of 11.2 ml. In this example, 5.2 ml would be added to
the original 6.0 ml for a final volume 11.2 ml. The required volume of
PBS-G was added to the (ICS) to make the Final Cell Suspension (FCS).
The LY# count and PBS-G volume were recorded on the SPECTROX™ Test
20 Worksheet.

After redilution, a 200 μ l aliquot of the rediluted cell suspension was transferred into a new 12 x 75 mm tube and a cell count performed as described. The rediluted cell suspension printout (Final Cell Suspension LY#) was attached to the Worksheet and the inoculation
5 volume recorded.

EXAMPLE 3

10 Assessment of Glutathione Concentration

Plate Inoculation

The final cell suspension was placed into a sterile trough and a microtiter plate containing media was placed inside a laminar flow hood. Using a 12-channel manual micropipettor equipped with sterile 0-
15 50 μ l barrier tips, the specified amount (according to TABLE I) of final cell suspension was dispensed to each well of the plate.

TABLE I

The following volumes are used for plate inoculation based on the listed final cell suspension Lymphocyte number (LY#).

Final Cell Suspension LY#	Adjusted Volume For Inoculation
3.9 - .37	8.0 μ l
3.6	8.5 μ l
2.5 - 3.5	10.0 μ l
2.4	12.5 μ l
2.3	13.0 μ l
2.2	14.0 μ l
2.1	14.5 μ l
2.0	15.0 μ l
1.9	16.0 μ l
1.8	17.0 μ l
1.7	18.0 μ l
1.6	19.0 μ l
1.5	20.0 μ l
1.4	21.5 μ l
1.3	23.0 μ l
≤ 1.2	25.0 μ l

5

After addition of the cells to the media, the plates were covered and placed in a CO₂ incubator maintained at 37° for 96 hours.

Labeling

10

All labeling procedures were performed in the Radioisotope Room. The tritiated thymidine (H³-TdR) working solution was removed

from the refrigerator and warmed to 37°C in a water bath. After 96 hours, the microtiter plates were removed from the incubator. The H³-TdR working solution was placed in a sterile trough and a 12-channel manual micropipettor equipped with 0-50 µl sterile barrier tips was used to dispense 10 µl of the H³-TdR working solution into each well of the microtiter plate. The plate was returned to the 37°C incubator for 24 hours. The date and initials of the technician performing the labeling was recorded on the sample log sheet.

10 Harvesting

All harvesting procedures were performed in the Radioisotope Room. A single glass fiber filter mat (Packard Part # 6005416) was labeled with the sample number using a #2 pencil. The vacuum pump was turned on and the drying oven set at 100°C. The distilled water carboy attached to the harvester was filled. The microtiter plates were removed from the incubator 24 hours after the addition of H³-TdR. The date and initials of the technician performing harvesting was recorded on the sample log sheet.

With the Cell Harvester (Packard Model # C9619) in the open position, O-Rings exposed, the glass fiber filter mat was placed onto the harvester with the rough side touching the O-Rings. The cell harvester

was closed and the filter mat was wet with distilled water from the rinse tray. The harvester was left on vacuum cycle (VAC). The lid was removed from microtiter plate and the plate was placed under harvester probe tips. The plate was slowly raised onto harvester probes, until the tips of the probe touched the bottom of the plate. With the media aspirated, the bottoms of the wells were scrubbed with the probe tips by moving the microtiter plates slowly in a circular motion. Scrubbing was continued for 10 seconds. With the microtiter plate in contact with the harvester probe tips, scrubbing was continued and the "WASH" button was pressed for 10 seconds. Liquid was aspirated from wells and the above scrubbing steps were repeated. The plate was removed, the rinse tray filled with methanol, the tray raised, the methanol aspirated and the tray lowered.

The harvester was opened, with the filter mat adhering to the upper section and continued to operate on VAC for 5 seconds. After 5 seconds, the VAC was turned off and filter mat was removed from the harvester surface simultaneously. The filter mat was placed rough side up in drying oven for 10 minutes. The filter mats were removed from the oven and cooled to room temperature.

20

Counting of Radioactivity

All counting procedures were conducted in the Radioisotope Room. The filter mat loaded into counting cassette, rough side up. A collimator, thin stainless steel plate that holds filter mat in place, was placed over the filter mat. The cassettes were loaded into Packard Matrix 9600 Beta Particle Radioactivity Counter, and the flow of Q-Gas (1.3% n-butane in helium) started into Matrix 9600. The "START" button was pressed activating the counting protocol, counting the total radioactivity in each well for 3 minutes. Each sample count was stored in the hard drive of the Matrix 9600. In addition, a hard copy of the raw radioactivity counts was printed out.

Data Transformation

Data was downloaded from the Matrix 9600 hard drive onto a 3.5 diskette. The raw data was transformed into a reportable format using a macro executed in Microsoft Excel. This macro subtracts the plate background from each data point, generates an average for the triplicate well values, and presents this value as a percentage of the Plate Control value which is set equal to 100%.

Data Analysis (Normalization)

The antioxidant function of glutathione (GSH) was examined by lymphocyte growth response assay. The cells were placed in a chemically-defined, serum-free medium which contained 5 μ M of L-
5 Buthionine-(S,R)-sulfoximine (BSO), and the cells were stimulated to grow by the addition of a mitogenic compound. Growth response, measured by ³H-thymidine incorporation, is expressed as percentage of control growth. BSO suppresses the synthesis of glutathione transferase which is required for glutathione synthesis. Thus, the addition of BSO to
10 the media and its suppression of glutathione synthesis is reflected by a change in cell growth, indicative of the status of glutathione in the lymphocyte. The doses of BSO used (5 μ M) were established in 1995.

The individual's lymphocyte growth is expressed as a percentage of growth in media with or without BSO. Microsoft Excel was
15 used to calculate the ratio of +BSO% growth to the 100% growth (without BSO). Over 1000 patient's result ratios are entered and statistical analysis was performed by standard Excel programs to determine mean, median, range and variance of the ratio. "Outliers", data that falls outside of a ± 2 standard deviation range, were eliminated
20 from the population base. The remaining data was used to determine the reference point and establish the population reference range. This

population reference range was represented graphically as scatter plot distribution. The normal/reference value for glutathione has been determined to be greater than 85% (>85 %); therefore, any person with a glutathione value of less than 85 % will be classified as having a glutathione deficiency.

Equipment and Reagents

The following equipment was used in the assay of the present invention: Laminar Flow Hood; Centrifuge, Beckman GS-6; Cell Counter (Coulter Model T540); 12 Channel Pipettor (5 - 500 μ l); Electric Pipet Pump (Drummond); Sterile 50 ml conical plastic tubes with caps; 12 x 75 mm polypropylene tubes; Sterile 15 ml conical plastic centrifuge tubes with caps; Pipettes (0-20, 0-200, 0-100 μ l ranges); Sterile Glass Disposable Pipettes - 5.0 ml, 10.0 ml; Test Tube Racks and Aerosol Barrier Pipet Tips (0 - 50 μ l).

TABLE II shows the various reagents and their sources used in the method of the present invention:

TABLE II

20

Reagents

Adenine Hydrochloride	Sigma	A 8751
Antibiotic Solution (PSF)	GIBCO	15245-012
Arginine Hydrochloride	Sigma	A 5131

L-Buthionine-[S.R]-Sulfoximine	Sigma	B2515
d-Biotin	Sigma	B 4501
Calcium Chloride, Anhydrous	Sigma	C 4901
Choline Chloride	Sigma	C 1879
Cumene Hydroperoxide	Sigma	C 0524
Cyanocobalamin (Vit. B ₁₂)	Sigma	V 2876
Cystein Hydrochloride, Anhyd.	Sigma	C 1276
Disodium EDTA	Sigma	E 4884
Ferrous Sulfate Heptahydrate	Sigma	F 8633
Folinic Acid, Calcium Salt	Sigma	F 7878
Glucose	Sigma	G 5767
Glucose Solution (10%)	Sigma	G 3126
L-Glutamine	Sigma	G 3126
Glycine	Sigma	G 7126
HEPES, Free Acid	Sigma	H 3375
L-Histidine(HCL monohydrate	Sigma	H 8125
Histopaque (Ficoll/Diatrizoate)	Sigma	1077-1
Hydroxocobalamin(HCI (B ₁₂)	Sigma	H 7126
myo-Inositol	Sigma	I 5125
L-Isoleucine	Sigma	I 2752
L-Leucine	Sigma	L 8000
L-Lysine	Sigma	L 5626
Magnesium Sulfate, Anhydrous	Sigma	M 7506
Methanol, Absolute	VWR	VWR4300-7
L-Methionine	Sigma	M 9625
Niacinamide (Vitamin B ₃)	Sigma	N 3376
D-Pantothenate, Calcium	Sigma	P 0290
Phenol Red Solution, 0.5% PBS	Sigma	P 0290
L-Phenylalanine	Sigma	P 2126
PBS pH7.4	Sigma	P 3813
Phytohemagglutinin PHA-P	Sigma	L 8754
Potassium Phosphate, Dibasic	Sigma	P3786
Pyridoxine(HCI (Vitamin B ₆)	Sigma	P 9755
Riboflavin (Vitamin B ₂)	Sigma	R 4500
L-Serine	Sigma	S 4500
Sodium Chloride	Sigma	S 9625
Sodium Hydroxide, 5.0N Solution	VWR	RS 415 1
Sodium Pyruvate	Sigma	P 2256
Sodium Pyruvate Sol. (100mM)	Sigma	S 8636

Thiamin (Vitamin B ₁)	Sigma	T 4625
L-Threonine	Sigma	T 8625
Thymidine	Sigma	T 9250
[³ H] - Thymidine	ICN	24066
L-Tryptophan	Sigma	T 0254
L-Tyrosine	Sigma	T 3754
L-Valine	Sigma	V 0500

Solutions

All solutions were prepared using tissue culture grade
5 deionized water (tcd H₂O).

1. Phosphate Buffered Saline: To prepare a phosphate
buffered saline (PBS)+ 0.72% Glucose (PBS-G) solution, PBS was
prepared according to package instructions using tcd H₂O. Sufficient
10 10% glucose solution was added to achieve final glucose concentration
of 0.72%. The solution was sterilized by filtration and stored in a
refrigerator at (2-8°C).

2. Stock Medium: A concentrated (2X) stock media was
prepared in the following way: (1) 23.80 g HEPES; (2) 14.02 g Sodium
Chloride; (3) 1.05 g Dibasic Potassium Phosphate; (4) 0.241 g
15 Magnesium Sulfate; (5) 1.0 ml (10 µM) Adenine Hydrochloride; (6) 30.0
ml (100 mM) Sodium Pyruvate; (7) 0.5 ml 0.5% Phenol Red; (8) 5.0 ml
Antibiotic Mixture; (9) 8.0 ml 5N Sodium Hydroxide; (10) 20.0 ml
Fe/EDTA (1.0 mM FeSO₄/0.4 mM Na₂EDTA) were combined in tcd H₂O.
After all materials had been mixed thoroughly, pH was adjusted to 7.60
20 using 5N Sodium Hydroxide. A final volume of 1.0 L was achieved with

ted H₂O. The solution was sterilized by filtration through a 0.2 µM filter and stored in a refrigerator at 4°C. Stability of the solution is 4 weeks.

3. Basal Medium: Basal Media used for 100% plate control and the novel method of the present invention is prepared as follows: After filtration, sterile technique should be employed. The Basal Media was prepared under laminar flow hoods, and the ingredients, as presented in Table III, were mixed and brought to the final desired volume with tcd H₂O. The solution was sterilized by vacuum filtration into sterile bottles. The proper volume of PHA Stock Solution was added.

TABLE III

	FINAL VOLUME BASAL MEDIA (ML)					
STOCK SOLUTION	100 ml	250 ml	500 ml	1000 ml	1500 ml	2000 ml
2x Stock Media (ml)	50	125	250	500	750	1000
Thiamin (B ₁) Stock (μl)	10	25	50	100	150	200
Riboflavin (B ₂) Stock (μl)	10	25	50	100	150	200
Niacinamide (B ₃) Stock (μl)	10	25	50	100	150	200
Pyridoxine (B ₆) Stock (μl)	10	25	50	100	150	200
Vitamin B12 Stock (μl)	10	25	50	100	150	200
2nd Stock Folinic (μl)	10	25	50	100	150	200
Pantothenate Stock (μl)	10	25	50	100	150	200
Biotin Stock (μl)	10	25	50	100	150	200
Stock Glucose (ml)	0.72	1.8	3.6	7.2	10.8	14.4
Stock Chol/Ino (ml)	1.0	2.5	5.0	10.0	15.0	20.0
Stock All Aminos (ml)	1.0	2.5	5.0	10.0	15.0	20.0
Stock CaCl ₂ (ml)	0.5	1.25	2.5	5.0	7.5	10.0
PHA (ml)	0.2	0.5	1.0	2.0	3.0	4.0

4. L-Buthionine-[S.R.]-Sulfoximine (BS) Working Solution

5 Preparation (5μM/L): 0.2223g BSO (Sigma B 2515) is dissolved in 20ml PBS-G to make a 50mM BSO stock solution. The solution is then sterilized by vacuum filtration into a sterile bottle using aseptic technique. A second BSO stock solution (500μM) is prepared by adding

100 μ l of 50mM BSO into 900 μ l sterile PBS-G. The 5 μ M BSO working solution is prepared by adding 2.5ml of the second, 500 μ M BSO stock solution to 247.5ml of Basal media and mixing well. The stability of this solution is for 1 week.

5 5. Thymidine(ThY) Stock Solution: The cold thymidine (ThY) stock solution (1.33 mM ThY, 0.322 g/L) was prepared as follows: 0.161 g ThY (Sigma T 9250) was weighed and dissolved in tcd H₂O to final volume of 500 ml. The solution was sterilized by vacuum filtration into a sterile bottle using aseptic technique. The solution was aliquoted
10 into 50 ml centrifuge tubes. For short-term storage, the solution can be refrigerated (4°C) with stability for one month. For long-term storage, the solution can be frozen (-70°C) with stability for 6 months. The thymidine working solution should be used to dilute radioactive thymidine (H³TdR) for labeling of cells.

15 To prepare the thymidine working solution, ThY + ³H-TdR, the following chemicals were added to 300 ml sterile, de-gassed tcd H₂O: 1.15 ml of 1.33 mM ThY (Cold), 1.70 ml of ³H-TdR (ICN part #24066) with a specific activity of 300 μ Ci/mmol. Once made, the solution may be stored under refrigerated conditions (4°C) for 1 week.

20

EXAMPLE 4

Assessment of Cysteine Concentration

Plate Inoculation

- 5 The final cell suspension was placed into a sterile trough and a microtiter plate containing media was placed inside a laminar flow hood. Using a 12-channel manual micropipettor equipped with sterile 0-50 μ l barrier tips, the specified amount (according to TABLE IV) of final cell suspension was dispensed to each well of the plate.

TABLE IV

The following volumes are used for plate inoculation based on the listed final cell suspension Lymphocyte number (LY#).

5

Final Cell Suspension LY#	Adjusted Volume For Inoculation
3.9 - .37	8.0 μ l
3.6	8.5 μ l
2.5 - 3.5	10.0 μ l
2.4	12.5 μ l
2.3	13.0 μ l
2.2	14.0 μ l
2.1	14.5 μ l
2.0	15.0 μ l
1.9	16.0 μ l
1.8	17.0 μ l
1.7	18.0 μ l
1.6	19.0 μ l
1.5	20.0 μ l
1.4	21.5 μ l
1.3	23.0 μ l
≤ 1.2	25.0 μ l

After addition of the cells to the media, 10 μ l of N-Acetyl-L-Cysteine (150 mM as final concentration) solution was added to columns 1-9 In three rows of the plates.

10

In columns 1-3 of another six rows of the plate, 10 μ l of 200 μ M Cumene Hydroperoxide (CuOOH) was added. In columns 4-6 on six

rows of the plate, 10 μ l of 300 μ M Cumene Hydroperoxide (CuOOH) was added. In columns 7-9 on six rows of the plate, 10 μ l of 400 μ M Cumene Hydroperoxide (CuOOH) was added. After addition of the CuOOH to the plates, the plates were covered and placed in a CO₂ incubator maintained at 37° for 96 hours.

Labeling

All labeling procedures were performed in the Radioisotope Room. The tritiated thymidine (H³-TdR) working solution was removed from the refrigerator and warmed to 37°C in a water bath. After 96 hours, the microtiter plates were removed from the incubator. The H³-TdR working solution was placed in a sterile trough and a 12-channel manual micropipettor equipped with 0-50 μ l sterile barrier tips was used to dispense 10 μ l of the H³-TdR working solution into each well of the microtiter plate. The plate was returned to the 37°C incubator for 24 hours. The date and initials of the technician performing the labeling was recorded on the sample log sheet.

Harvesting

All harvesting procedures were performed in the Radioisotope Room. A single glass fiber filter mat (Packard Part #

6005416) was labeled with the sample number using a # 2 pencil. The vacuum pump was turned on and the drying oven set at 100°C. The distilled water carboy attached to the harvester was filled. The microtiter plates were removed from the incubator 24 hours after the
5 addition of H³-TdR. The date and initials of the technician performing harvesting was recorded on the sample log sheet.

With the Cell Harvester (Packard Model # C9619) in the open position, O-Rings exposed, the glass fiber filter mat was placed onto the harvester with the rough side touching the O-Rings. The cell harvester
10 was closed and the filter mat was wet with distilled water from the rinse tray. The harvester was left on vacuum cycle (VAC). The lid was removed from microtiter plate and the plate was placed under harvester probe tips. The plate was slowly raised onto harvester probes, until the tips of the probe touched the bottom of the plate. With the media
15 aspirated, the bottoms of the wells were scrubbed with the probe tips by moving the microtiter plates slowly in a circular motion. Scrubbing was continued for 10 seconds. With the microtiter plate in contact with the harvester probe tips, scrubbing was continued and the "WASH" button was pressed for 10 seconds. Liquid was aspirated from wells and the
20 above scrubbing steps were repeated. The plate was removed, the rinse

tray filled with methanol, the tray raised, the methanol aspirated and the tray lowered.

The harvester was opened, with the filter mat adhering to the upper section and continued to operate on VAC for 5 seconds. After 5 seconds, the VAC was turned off and filter mat was removed from the harvester surface simultaneously. The filter mat was placed rough side up in drying oven for 10 minutes. The filter mats were removed from the oven and cooled to room temperature.

10 Counting of Radioactivity

All counting procedures were conducted in the Radioisotope Room. The filter mat loaded into counting cassette, rough side up. A collimator, thin stainless steel plate that holds filter mat in place, was placed over the filter mat. The cassettes were loaded into Packard Matrix 9600 Beta Particle Radioactivity Counter, and the flow of Q-Gas (1.3% n-butane in helium) started into Matrix 9600. The "START" button was pressed activating the counting protocol, counting the total radioactivity in each well for 3 minutes. Each sample count was stored in the hard drive of the Matrix 9600. In addition, a hard copy of the raw radioactivity counts was printed.

Data Transformation

Data was downloaded from the Matrix 9600 hard drive onto a 3.5 diskette. The raw data was transformed into a reportable format using a macro executed in Microsoft Excel. This macro subtracts the plate background from each data point, generates an average for the triplicate well values, and presents this value as a percentage of the Plate Control value which is set equal to 100%.

Data Analysis (Normalization)

This test is an assessment of intracellular cysteine concentration, which is one determinant of a cell's antioxidant capacity. Using lymphocytes stimulated to grow by mitogen, the antioxidant function is expressed by comparing the difference in growth response of lymphocytes, with or without cysteine, in the presence of Cumene Hydroperoxide (CuOOH). The CuOOH creates the oxidative stress used to measure overall antioxidant status in individual patient. The addition of cysteine to the media provides an antioxidant with the ability to repair damage caused by the oxidative stress caused by CuOOH. The doses of CuOOH and cysteine were 300 μ M and 150 μ M, respectively.

The individual's lymphocyte growth expressed as a percentage of growth in media with CuOOH or cysteine + CuOOH.

Microsoft Excel was used to calculate the ratio of Cysteine + CuOOH growth to the CuOOH % growth. More than 1000 patient's result ratios were entered for statistical analysis by standard Excel programs to determine the mean, median, range and variance. Outliers--data that
5 falls outside of a ± 2 standard deviation range--were eliminated from the population base. The remaining data was used to determine the reference point and establish the population reference range. This population reference range was graphically represented as scatter plot distribution. The normal/reference value for cysteine has been
10 determined to be less than 127% (<127%), therefore, any person with a cysteine value of 127% or greater will be classified as having a cysteine deficiency.

Equipment and Reagents

15 The following equipment was used in the assay of the present invention: Laminar Flow Hood; Centrifuge, Beckman GS-6; Cell Counter (Coulter Model T540); 12 Channel Pipettor (5 - 500 μ l); Electric Pipet Pump (Drummond); Sterile 50 ml conical plastic tubes with caps; 12 x 75 mm polypropylene tubes; Sterile 15 ml conical plastic centrifuge
20 tubes with caps; Pipettes (0-20, 0-200, 0-100 μ l ranges); Sterile Glass

Disposable Pipettes - 5.0 ml, 10.0 ml; Test Tube Racks and Aerosol Barrier Pipet Tips (0 - 50 μ l).

TABLE V shows the various reagents and their sources used in the method of the present invention:

5

TABLE V

Reagents

Adenine Hydrochloride	Sigma	A 8751
Antibiotic Solution (PSF)	GIBCO	15245-012
Arginine Hydrochloride	Sigma	A 5131
N-Acetyl-L-Cysteine	Sigma	A8199
d-Biotin	Sigma	B 4501
Calcium Chloride, Anhydrous	Sigma	C 4901
Choline Chloride	Sigma	C 1879
Cumene Hydroperoxide	Sigma	C 0524
Cyanocobalamin (Vit. B ₁₂)	Sigma	V 2876
Cystein Hydrochloride, Anhyd.	Sigma	C 1276
Disodium EDTA	Sigma	E 4884
Ferrous Sulfate Heptahydrate	Sigma	F 8633
Folinic Acid, Calcium Salt	Sigma	F 7878
Glucose	Sigma	G 5767
Glucose Solution (10%)	Sigma	G 3126
L-Glutamine	Sigma	G 3126
Glycine	Sigma	G 7126
HEPES, Free Acid	Sigma	H 3375
L-Histidine(HCL monohydrate	Sigma	H 8125
Histopaque (Ficoll/Diatrizoate)	Sigma	1077-1
Hydroxocobalamin(HCl (B ₁₂)	Sigma	H 7126
myo-Inositol	Sigma	I 5125
L-Isoleucine	Sigma	I 2752
L-Leucine	Sigma	L 8000
L-Lysine	Sigma	L 5626
Magnesium Sulfate, Anhydrous	Sigma	M 7506
Methanol, Absolute	VWR	VWR4300-7
L-Methionine	Sigma	M 9625
Niacinamide (Vitamin B ₃)	Sigma	N 3376

D-Pantothenate, Calcium	Sigma	P 0290
Phenol Red Solution, 0.5% PBS	Sigma	P 0290
L-Phenylalanine	Sigma	P 2126
PBS pH7.4	Sigma	P 3813
Phytohemagglutinin PHA-P	Sigma	L 8754
Potassium Phosphate, Dibasic	Sigma	P3786
Pyridoxine(HCl (Vitamin B ₆))	Sigma	P 9755
Riboflavin (Vitamin B ₂)	Sigma	R 4500
L-Serine	Sigma	S 4500
Sodium Chloride	Sigma	S 9625
Sodium Hydroxide, 5.0N Solution	VWR	RS 415 1
Sodium Pyruvate	Sigma	P 2256
Sodium Pyruvate Sol. (100mM)	Sigma	S 8636
Thiamin (Vitamin B ₁)	Sigma	T 4625
L-Threonine	Sigma	T 8625
Thymidine	Sigma	T 9250
[³ H] - Thymidine	ICN	24066
L-Tryptophan	Sigma	T 0254
L-Tyrosine	Sigma	T 3754
L-Valine	Sigma	V 0500

Solutions

All solutions were prepared using tissue culture grade deionized water (tcd H₂O).

- 5 1. Phosphate Buffered Saline: To prepare a phosphate buffered saline (PBS)+ 0.72% Glucose (PBS-G) solution, PBS was prepared according to package instructions using tcd H₂O. Sufficient 10% glucose solution was added to achieve final glucose concentration of 0.72%. The solution was sterilized by filtration and stored in a 10 refrigerator at (2-8°C).

2. Stock Medium: A concentrated (2X) stock media was prepared in the following way: (1) 23.80 g HEPES; (2) 14.02 g Sodium Chloride; (3) 1.05 g Dibasic Potassium Phosphate; (4) 0.241 g Magnesium Sulfate; (5) 1.0 ml (10 μ M) Adenine Hydrochloride; (6) 30.0 ml (100 mM) Sodium Pyruvate; (7) 0.5 ml 0.5% Phenol Red; (8) 5.0 ml Antibiotic Mixture; (9) 8.0 ml 5N Sodium Hydroxide; (10) 20.0 ml Fe/EDTA (1.0 mM FeSO_4 /0.4 mM Na_2EDTA) were combined in tcd H_2O . After all materials had been mixed thoroughly, pH was adjusted to 7.60 using 5N Sodium Hydroxide. A final volume of 1.0 L was achieved with tcd H_2O . The solution was sterilized by filtration through a 0.2 μ M filter and stored in a refrigerator at 4°C. Stability of the solution is about 4 weeks.

3. Basal Medium: Basal Media used for 100% plate control and the novel method of the present invention is prepared as follows: After filtration, sterile technique should be employed. The Basal Media was prepared under laminar flow hoods, and the ingredients, as presented in Table VI, were mixed and brought to the final desired volume with tcd H_2O . The solution was sterilized by vacuum filtration into sterile bottles. The proper volume of PHA Stock Solution was added.

TABLE VI

	FINAL VOLUME BASAL MEDIA (ML)					
STOCK SOLUTION	100 ml	250 ml	500 ml	1000 ml	1500 ml	2000 ml
2x Stock Media (ml)	50	125	250	500	750	1000
Thiamin (B ₁) Stock (μl)	10	25	50	100	150	200
Riboflavin (B ₂) Stock (μl)	10	25	50	100	150	200
Niacinamide (B ₃) Stock (μl)	10	25	50	100	150	200
Pyridoxine (B ₆) Stock (μl)	10	25	50	100	150	200
Vitamin B12 Stock (μl)	10	25	50	100	150	200
2nd Stock Folinic (μl)	10	25	50	100	150	200
Pantothenate Stock (μl)	10	25	50	100	150	200
Biotin Stock (μl)	10	25	50	100	150	200
Stock Glucose (ml)	0.72	1.8	3.6	7.2	10.8	14.4
Stock Chol/Ino (ml)	1.0	2.5	5.0	10.0	15.0	20.0
Stock All Aminos (ml)	1.0	2.5	5.0	10.0	15.0	20.0
Stock CaCl ₂ (ml)	0.5	1.25	2.5	5.0	7.5	10.0
PHA (ml)	0.2	0.5	1.0	2.0	3.0	4.0

4. Cumene Hydroperoxide (CuOOH) Stock Solutions;

5 Cumene hydroperoxide (Sigma C 0524) has a limited shelf life. Expiration date of the material is three (3) months from date of receipt from Sigma. When pipetting CuOOH from the bottle, one should be aware of the increased viscosity of CuOOH. One must ensure that

aspiration of liquid into the pipet tip is complete. This requires extra time and attention to ensure adequate fill of pipet tip. Also, one must wipe the tip to remove excess CuOOH on the outside of the pipet tip. Likewise, one must dispense the CuOOH completely. Following this
5 procedure, one may store the solution under refrigerated conditions (2-8°C) with stability for up to 3 months.

For the first CuOOH stock solution (1.0 M CuOOH in PBS-G), 9.5 µl cumene hydroperoxide (CuOOH) are mixed with 990.5 µl PBS-G. Following this procedure, one may store the solution under refrigerated
10 conditions (2-8°C) with stability for up to 3 months.

The second CuOOH stock solution (100 mM in PBS-G) must be prepared daily before cell isolation, and should not be stored overnight. To prepare the second CuOOH stock solution, 200 µl of the first stock solution were mixed with 1800 µl of PBS-G.

15 For the cumene hydroperoxide working solutions, 4 working solutions were prepared daily before cell isolation. The solution was added to the CuOOH Transfer Plate (separate microtiter plate) for loading of patient plates. The solutions should not be stored overnight. The working solution was 300 µM CuOOH: 330 µl of the second CuOOH
20 stock solution was mixed with 4670 µl PBS-G.

5. Thymidine(ThY) Stock Solution: The cold thymidine (ThY) stock solution (1.33 mM ThY, 0.322 g/L) was prepared as follows: 0.161 g ThY (Sigma T 9250) was weighed and dissolved in tcd H₂O to final volume of 500 ml. The solution was sterilized by vacuum filtration
5 into a sterile bottle using aseptic technique. The solution was aliquoted into 50 ml centrifuge tubes. For short-term storage, the solution can be refrigerated (4°C) with stability for one month. For long-term storage, the solution can be frozen (-70°C) with stability for 6 months. The thymidine working solution should be used to dilute radioactive
10 thymidine (H³TdR) for labeling of cells.

To prepare the thymidine working solution, ThY + ³H-TdR, the following chemicals were added to 300 ml sterile, de-gassed tcd H₂O: 1.15 ml of 1.33 mM ThY (Cold), 1.70 ml of ³H-TdR (ICN part #24066), and 300 µCi/mmol specific activity. Once made, the solution may be
15 stored under refrigerated conditions (4°C) for 1 week.

6. N-Acetyl-L-Cysteine (NAC) Working Solution: 0.0702g of N-Acetyl-L-Cysteine (Sigma A8199) was dissolved in 4ml PBS-G solution to form a final concentration 0.1M stock solution. This stock
20 solution was dispensed in 1 ml aliquots into 12 x 75 sterile tubes, and

these aliquots may be stored up to one week in the freezer at -70°C. To make the working solution, take 990µl of one aliquot of the 0.1M NAC stock solution and add to 29.01ml of the PBS-G solution and mix well. The working solution should be sterilized by vacuum filtration into a sterile bottle, using aseptic technique. Finally, 10µl of this solution should be added into the triplicate test wells. It is important that the N-Acetyl-L-Cysteine working solution should be made fresh every day before the test is performed.

One skilled in the art will appreciate readily that the present invention is well adapted to carry out the objects and obtain the ends and advantages mentioned, as well as those objects, ends and advantages inherent herein. The present examples, along with the methods, procedures, treatments, molecules, and specific compounds described herein are presently representative of preferred embodiments, are exemplary, and are not intended as limitations on the scope of the invention. Changes therein and other uses will occur to those skilled in the art which are encompassed within the spirit of the invention as defined by the scope of the claims.